





Danis Seal

Esophageal Stent Danis Seal

Removable fully covered nitinol stent

Sealing esophageal leaks and ruptures

Large diameter with optimum length
and expansion force



Designed & manufactured
in the **Czech Republic (EU)**

www.ellacs.eu

Manufacturer:
ELLA-CS, s.r.o.
Milady Horakove 504/45, Trebes
500 06 Hradec Kralove, Czech Republic

Indication

Treatment of esophageal or gastric leakage from:


- Anastomotic wound dehiscence (e.g. esophagectomy, gastric bypass)
- Esophageal ruptures / perforations (spontaneous rupture - Boerhaave's syndrome; iatrogenic rupture / perforation occurring during esophageal dilations, endoscopic manipulations, traumatic esophageal ruptures due to blunt thoracic trauma)

Features and Benefits

- + **Variable stent braiding** – reduced migration.
- + **Radiopaque markers at both stent ends and at the midpoint** – accurate stent positioning.
- + **Metallic retrieval loops** – simple endoscopic stent removal with ELLA Extractor.
- + **Durable silicone covering** – high mechanical durability and acid resistance.

Available sizes

REF No.	Esophageal Stent Danis Seal			Delivery system	
	Stent flares diameter [mm]	Stent body diameter [mm]	Nominal length [mm]	Active length [cm]	Outer diameter [F]
019-11-25-135	30	25	135	75	28 / 18



PUBLICATIONS



www.ellacs.eu
info@ellacs.eu

CE 1014

Manufacturer:
ELLA-CS, s.r.o.
Milady Horakove 504/45, Trebes
500 06 Hradec Kralove, Czech Republic

SX-11-AJ-08/13/P/REV-3-04/2024