




Danis Stent

Procedure Pack with SX-ELLA Stent Danis

Worldwide unique solution of esophageal variceal bleeding

Removable fully covered nitinol stent

Complete procedure pack for emergency use



Designed & manufactured
in the Czech Republic (EU)

www.ellacs.eu

3D ANIMATION



Manufacturer:
ELLA-CS, s.r.o.
Milady Horakove 504/45, Trebes
500 06 Hradec Kralove, Czech Republic

Indication

SX-ELLA Stent Danis is intended to stop acute and/or refractory bleeding from esophageal varices.

Additional description

Danis Procedure Pack (019-08S-25-135) consists of a plastic carry case containing all items required for stent implantation:

- Delivery system preloaded with stent, guide wire, plastic syringes 50 ml and 20 ml, bite block, paper emesis basin, slobber cloth, waste bag with plastic clip.

Danis Procedure Pack - Basic (019-08S-25-135-B) consists of:

- Delivery system preloaded with stent, guide wire, plastic syringe 50 ml.

The delivery system 9.4 mm (28 F) / 6.6 mm (20 F) is available in a standard length of 60 cm.

Available sizes


REF No.	SX-ELLA Stent Danis			Delivery system	
	Stent flares diameter [mm]	Stent body diameter [mm]	Nominal length [mm]	Active length [cm]	Outer diameter [F]
019-08S-25-135	30	25	135	60	28 / 20
019-08S-25-135-B	30	25	135	60	28 / 20

Features and Benefits

- + **Unique delivery system** – possible implantation without endoscopy / X-ray control.
- + **Early discharge from ICU and open stent lumen** – improved patient quality of life compared to balloon tamponade¹.
- + **Variable stent braiding** – reduced migration rates.
- + **Metallic retrieval loops** – easy endoscopic stent removal within 7 days.
- + **Significantly reduced complications and SAEs** when compared to balloon tamponade in acute oesophageal varices².

¹ NICE Medical technologies guidance on www.nice.org.uk/guidance/mtg57

² Escorsell A, et al., 'Esophageal Balloon Tamponade Versus Esophageal Stent in Controlling Acute Refractory Variceal Bleeding: A Multicenter Randomized Controlled Trial' (2016) Hepatology, (Vol 63, No.6)



PUBLICATIONS



www.ellacs.eu
info@ellacs.eu

CE 2797

Manufacturer:
ELLA-CS, s.r.o.
Milady Horakove 504/45, Trebes
500 06 Hradec Kralove, Czech Republic

SX-08-AJ-11/05/P/REV-12-04/2024